

# ScanX<sup>®</sup> Duo

DIGITAL IMAGING SYSTEM

## ANNOUNCING SCANX<sup>®</sup> DUO

THE NEW COMPACT INTRAORAL SCANX PACKED WITH PERFORMANCE ADVANTAGES.

**ScanX Duo works the way you do.** You want your digital system to be simple to operate and quick to master. Duo builds on the imaging techniques already used by your staff and minimizes adjustment of routines. You won't even need new X-ray equipment.

**Intelligent Track Control for chairside image processing.** An imaging first and ALLPRO exclusive advantage! Each of the two processing slots can be simultaneously and independently controlled. One Duo can be used by two operators at the same time with different patients.

**In-Line Erase.** Duo will read the imaging PSP (Phosphor Storage Plate) and then automatically erase it in the same continuous cycle. In-Line Erase saves time and simplifies the imaging process.

**Duo uses a superior receptor, a Phosphor Storage Plate.** Imagine digital film that is ultra-thin, wireless, flexible and reusable thousands of times. Plus it is very inexpensive. Plates allow accurate positioning within the oral cavity, minimizing the chance of image elongation, foreshortening or retakes. That's superior technology.

Understand your choices so you can make an informed decision. Arrange for a complete ScanX demonstration at your practice. Learn why ScanX is the most versatile digital imaging technology available.



- ▶ Intelligent Track Control (ITC) enables independent, yet simultaneous control of each processing slot.
- ▶ Duo is so small, just 11 3/4" H x 10" W x 10 1/2" D.

 **ALLPRO** IMAGING  
THE IMAGING EXPERTS

# scanX duo

DIGITAL IMAGING SYSTEM



## SPECIFICATIONS

Dimensions:	11 3/4" H x 10" W x 10 1/2" D
Weight:	25 lbs.
Electrical:	100 - 240V AC, 50/60 Hz
Resolution:	Up to 18 lp/mm (depending on scanning model/pixel size selected)
Ports:	USB
Laser Classification:	Class I Laser Product, Compliance with FDA 21 CFR 1040.10 and IEC 60825-1
Imaging Medium:	Two processing slots accept Phosphor Storage Plates. Intraoral sizes 0, 1, 2, 3, 4
Accessories Included:	Power Cord, Computer Connector (USB), 2 each size 2 PSP Guides, 1 each size 0 and 1 PSP Guides, 1 PSP Table Mat, Drivers, Utilities and Installation Manual CD
Warranty:	ScanX Duo is warranted to be free from defects for a period of 24 months.

## MINIMUM COMPUTER REQUIREMENTS / SPECIFICATIONS

CPU Speed:	2.8 GHz Pentium IV
Operating System:	Windows 2000 Professional with Service Pack 4 or later, or Windows XP Professional with Service Pack 1 and the KB B22603 update or later.
RAM:	1 GB
Hard Drive:	120 GB, 200 MB free disk space required to start screening
Monitor:	SVGA 17", 1024 x 768 or higher resolution, contrast ratio 450:1, .22 dot pitch capability. For optimum viewing a CRT is recommended.
Video Display Adapter:	32 MB RAM
Keyboard / Mouse:	Standard
Backup Device:	Recommended for daily archiving
External Surge Protector:	Recommended
Printer:	Recommended
USB Port:	Must be 1.1, 2.0 or later



For more information or to contact an authorized ALLPRO dealer near you, call 1-800-247-8326, Ext 5611 or visit [www.allproimaging.com](http://www.allproimaging.com)



1295 Walt Whitman Road, Melville, New York 11747-3062  
516-433-7676 • Fax: 516-433-7683 • [www.allproimaging.com](http://www.allproimaging.com)