

# Flat Top Surgery Table

By Technidyne



## Flat Top Surgery Table With Electric Height Adjustment

- Raises and lowers electrically with foot pedal
- Tilts 0 to 30° with fingertip lever control
- Uses rechargeable battery power to eliminate cords on floor
- Designed to be convenient for sit-down surgery and stand-up surgery
- 4 quick release restraints are attached to underside of table and can be easily positioned and locked
- Heater with digital temperature readout *available*, ( Battery operation not available with heater option)
- Drain hole and bucket *available*
- 3 year warranty

**TECHNIDYNE CORPORATION**

[www.technivet.com](http://www.technivet.com)

**1-800-654-8073**